

CARDIO PRESTIGE II STETHOSCOPE DUAL



Ref. 506011



13/06/2025 22:36:27 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

CARDIO PRESTIGE II STETHOSCOPE DUAL



Ref. 506011

With its heavy horn with "premium acoustic" technology, the new Cardio Prestige II is the most efficient stethoscope.

Description

Developed with the help of the CNRS and the CENTRALE school, the CARDIOPRESTIGE II stethoscope is the quintessence of SPENGLER know-how. With its heavy roof equipped with "premium acoustic" technology, its oversize tubing, it offers a level of performance never reached.

EXCEPTIONAL

- ACOUSTICS• Large diameter lyres•
- Multi-frequency floating membrane- Light pressure = listening to low frequencies (vascular and cardiac)
- Firm pressure = listening to high frequencies (pulmonary)

COMFORT•

- Anti-cold ring•
- Nano silicone ear tips•
- New lyre design more adapted to anatomy•
- Pavilion indicator

VERSATILITY

- Adult and pediatric pavilion

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446034399	HS Code	90181990
Warranty (months)	84		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.601
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

📅 13/06/2025 22:36:27 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

🌐 www.gsh-med.fr

☎ +33(0)4 42 90 31 31

✉ info@gsh-med.fr