

## CARDIO PRESTIGE II STETHOSCOPE SINGLE



Ref. 506010



14/06/2025 00:25:52 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## CARDIO PRESTIGE II STETHOSCOPE SINGLE



Ref. 506010

**With its heavy horn with "premium acoustic" technology, the Cardio Prestige II is the most efficient stethoscope.**

### Description

Developed with the help of the CNRS and the CENTRALE school, the CARDIOPRESTIGE II stethoscope is the quintessence of SPENGLER know-how. With its heavy roof equipped with "premium acoustic" technology, its oversize tubing, it offers a level of performance never reached.

#### EXCEPTIONAL

ACOUSTICS • Large diameter lyres •

Multi-frequency floating membrane- Light pressure = listening to low frequencies (vascular and cardiac)

-Firm pressure = listening to high frequencies (pulmonary)

COMFORT • Anti-cold ring •

Nano silicone ear tips

• New lyre design more adapted to anatomy

Available in double pavilion (adult and pediatric pavilion).

### Features

CE Marking	CE	Medical device	Class I
Bar code	3700446029975	HS Code	90181990
Warranty (months)	84		
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.543
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

📅 14/06/2025 00:25:52 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

🌐 [www.gsh-med.fr](http://www.gsh-med.fr)

☎ +33(0)4 42 90 31 31

✉ [info@gsh-med.fr](mailto:info@gsh-med.fr)