



MAGISTER STETHOSCOPE FOR ADULT / BLACK



Ref. 507200



09/01/2025 23:06:41 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

MAGISTER STETHOSCOPE FOR ADULT / BLACK



Ref. 507200

=====

Description

For general-purpose auscultation. Stainless steel dual head chestpiece for adults. Multifrequency diaphragm. Angled binaural with built-in spring and flexible, screw-on tips. Ideal grip and high performance make it the perfect product for general and specialist use.

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446009069	HS Code	90181990
Warranty (months)	60		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	16	Packed unit / Height (cm)	4
Packed unit / Length (cm)	32	Packed unit / Weight	0.398
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Colors	Black	Packed unit / Quantity	1

09/01/2025 23:06:41 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr