

DUAL PULSE STETHOSCOPE / BLACK

Ref. 507430



14/06/2025 00:06:30 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

DUAL PULSE STETHOSCOPE / BLACK



Ref. 507430

=====

Description

Pearly tube model. For daily blood pressure monitoring. Powder-coated or chromium-plated (black model) zinc alloy single head chestpiece. Epoxy resin, high resolution acoustic membrane. Screw-on flexible tips to match the PVC tubing. Angled binaural with built-in spring.

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446026530	HS Code	90181990
Warranty (months)	36		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.237
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Colors	Red	Packed unit / Quantity	1

14/06/2025 00:06:30 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr