

TEMPO EAR PRO + - EAR THERMOMETER (RECHARGEABLE)



Ref. THE003



10/07/2025 20:25:48 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

TEMPO EAR PRO + - EAR THERMOMETER (RECHARGEABLE)



Ref. THE003

Performance, hygiene and comfort for healthcare professionals

Description

Spengler's Tempo Ear Pro is the infrared ear thermometer designed for demanding medical environments. Fast, reliable and hygienic, it makes daily temperature taking easier, whether in the office, hospital or clinic.

- Ultra-fast measurement: less than 3 seconds
- Probe preheated to 36°C: patient comfort
- Infrared technology: reliable, even when moving
- Probe cover ejection: maximum hygiene without contact
- Memory of last 30 time-stamped measurements
- Automatic sleep mode after 60 seconds
- Pulse counter
- Lanyard
- Rechargeable

Included with Tempo Ear Pro + :

- 1 Tempo Ear Pro + ear thermometer
- 1 storage and charging cradle
- 20 probe covers
- 1 attachment cord
- 1 user manual
- Thermometer calibration and verification certificate with batch number

Also available : Tempo EAR PRO with batteries instead of battery.

Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	6970392219448	HS Code	9025199010
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	5	Packed unit / Height (cm)	3
Packed unit / Length (cm)	14	Packed unit / Weight	nc
Bare Unit / Weight (kg)	0.117		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	50

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr